

# AFFINITY PURIFICATION OF SOLUBLE GST-TAGGED PROTEINS

[www.ptglab.com](http://www.ptglab.com)

---

- 1. Lyse cells:**
  - a. Suspend the cell pellet in 30–35 ml of glutathione S-transferase (GST)-washing buffer with 10 mM PMSF and 0.5 M EDTA.
  - b. Sonicate cells in an ice-bath at 200 W for 6 min.
  - c. Rotate the lysed solution for 1 h at 4°C.
  - d. Centrifuge the cell lysate for approximately 13 min at 8000 rpm, 4°C.

---

- 2. Bind protein to beads:**
  - a. Transfer the supernatant to 600 µl of GST-beads.
  - b. Rotate the mixture overnight at 4°C.
  - c. Collect the beads by centrifugation at 2000 rpm for 10–30 seconds, 4°C. Collect the protein bound beads in eppendorf tubes.

---

- 3. Wash out the unbound proteins from beads:**
  - a. Wash the beads 3 times with 1ml of GST-washing buffer. Discard the supernatant.

---

- 4. Elute proteins from beads:**
  - a. Add 300 µl of GST-elution buffer to the beads.
  - b. Rotate the mixture for 1 h at 4°C.
  - c. Collect the supernatant by centrifugation at 300 rpm for 10–30 seconds.
  - d. Repeat steps 4 a–c.
  - e. Combine the eluent (total volume of 600 µl).
  - f. Check the molecular weight and purity of the enriched protein by SDS-PAGE analysis.

# AFFINITY PURIFICATION OF SOLUBLE GST-TAGGED PROTEINS

[www.ptglab.com](http://www.ptglab.com)

---

## Buffers Needed

<b>GST-washing buffer (PBST buffer)</b>	<b>1000 ml</b>
58 mM Na <sub>2</sub> HPO <sub>4</sub>	8.24 g
17 mM NaH <sub>2</sub> PO <sub>4</sub>	2.04 g
68 mM NaCl	3.98 g
1% Triton X-100	10 ml
Add ddH <sub>2</sub> O to 1000 ml	
Adjust to pH 7.4	

  

<b>GST-elution buffer</b>	<b>1000 ml</b>
100 mM GSH	30.70 g
10% Glycerol	100 ml
1x PBST buffer	900 ml
Adjust to pH 8.0	