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HUMANKINE® Protein Reconstitution Buffer

Types of Protein Reconstitution Buffer

1. Sterile water
2. Sterile 1 × PBS
3. Sterile 4 mM HCl
4. Sterile 10 mM acetic acid
5. Sterile 20% ethanol + 50 mM sodium acetate + 75 mM acetic acid

Note: Select the Reconstitution buffer based on the datasheet of each protein.

Product Description

1. 1 × PBS phosphate buffer (sterile) formulation: 154 mM NaCl, 1mM KH_2PO_4 , 5.58mM Na_2HPO_4 aqueous solution.

Storage: pH 7.2-7.4, stable for one month at room temperature, three months at 4 °C , two years at -20°C .

If crystals have formed at 4 °C or -20 °C , please balance to room temperature, stir until all the crystals are completely dissolved before use.

2. Sterile 4 mM HCl
Sterile 10 mM acetic acid
Sterile 20% ethanol + 50 mM sodium acetate + 75 mM acetic acid

Note: All these three protein reconstitution buffers **should be prepared freshly** due to containing of volatile components, long-term storage is not recommended.

Notice

1. This product is sterile and can be used in aseptic experiments. This product is preservative-free, please avoid microbial contamination, screw the cap tight immediately after use. Do not use if turbid insoluble component appears during storage.
2. 1 × PBS phosphate buffer (sterile) contains a lower concentration of phosphate. Bivalent or trivalent metal ions should be avoided during use, otherwise phosphate precipitates may form.
3. Do not vortex. Cytokines or recombinant proteins can be completely dissolved by gently pipetting several times. For insoluble or slow dissolving proteins, shake slowly on a shaker or standing the reconstitution solution for at least 2 hours at 4 °C . For insoluble cytokines, please refer datasheet.
4. **This product is free of protein protectant**, if long-term storage is needed, it is recommended to be stored at protective solution of 0.1% HSA or BSA (endotoxin-free).

For technical support and original validation data for this product please contact:

Proteintech Group, USA,
5400 Pearl Street, Suite 300,
Rosemont, IL 60018, USA
t. 1-888-478-4522
e. proteintech@ptglab.com

Proteintech Europe,
Manchester Science Park, Kilburn House,
Lloyd Street North, Manchester, M15 6SE
t. (+44)-161-22-66-144
e. europe@ptglab.com

Proteintech China,
D3-3, No.666 Gaoxin Avenue, Wuhan East Lake
Hi-tech Development Zone, Wuhan, P.R.C.
t. 86-27-87531629
e. Proteintech-CN@ptglab.com