

Protocol for removal of BSA from antibodies

BSA depletion agarose facilitates fast & easy removal of BSA from antibodies with high antibody recovery. BSA is often added to antibodies for stabilization but may interfere with conjugation methods and other downstream applications.

General Notes

- BSA Depletion Agarose is a high-affinity binding resin binding BSA to deplete it from an antibody solution. The antibody of interest will flow through.
- BSA Depletion Agarose is compatible with antibodies from a wide range of species such as rabbit, mouse, rat, chicken, sheep etc.
- BSA Depletion Agarose is also compatible with typical antibody buffer compounds, including glycerol up to 50%.
- BSA concentrations in antibodies can vary from 0.01% to 1%.

Our recommendation for usage:

- Use 50 μ l **BSA Depletion Agarose** (slurry) for 100 μ l of antibody supplemented with $\leq 0.1\%$ BSA.
Note: We do not recommend using less than 50 μ l to avoid significant loss of antibody due to bead handling issues.
- Use 200 μ l beads (slurry) for 100 μ l antibody solutions with 0.5% BSA.
- Use 300 μ l beads (slurry) for 100 μ l antibody solutions with 1% BSA.

Bead volumes can be scaled up proportionally to the antibody volume.

- Important note: Commercial BSA preparations often contain bovine serum proteins, including bovine immunoglobulins (IgG). As a result, the recovered IgG may appear lower in quantity, when in fact the purification has removed contaminating bovine IgG and yielded a cleaner preparation of the target IgG.

Procedure

Equilibration of beads (for 100 µl antibody volume)

1. Resuspend beads by pipetting up and down or by inverting. **Do not vortex.**
2. Remove stopper from a spin column and insert it into a 2 ml Eppendorf tube.
3. Transfer 100 µl (for BSA concentrations < 0.5%) or 300 µl (for BSA 1%) beads (slurry) to a spin column. Centrifuge at 2000x g for 1 min to remove storage solution.
4. Add 500 µl PBS to the beads and centrifuge at 2000x g for 1 min.
Discard the flow-through.

Depletion of BSA

5. Close the spin column and add your antibody (100 µl) to the equilibrated beads.
6. Incubate for 15 min at room temperature at 1000 rpm in a thermal shaker.
7. Open spin column and transfer to fresh 2 ml Eppendorf tube.
8. Centrifuge at 2000x g for 2 min to recover antibody.
9. Optional (may increase antibody yield): Add 100 µl PBS to the beads and centrifuge at 2000x g for 2 min. Combine eluate with eluate from step 8.

Note: This will decrease the concentration of the antibody.