

Protocol for removal of BSA from antibodies

BSA Depletion Magnetic Agarose facilitates fast & easy removal of BSA from antibodies with high antibody recovery. BSA is often added to antibodies for stabilization but may interfere with conjugation methods and other downstream applications.

General Notes

- BSA Depletion Magnetic Agarose is a high-affinity binding resin binding BSA to deplete it from an antibody solution. The antibody of interest will flow through.
- BSA Depletion Magnetic Agarose is compatible with antibodies from a wide range of species such as rabbit, mouse, rat, chicken, sheep etc.
- BSA Depletion Magnetic Agarose is also compatible with typical antibody buffer compounds, including glycerol up to 50%.
- BSA concentrations in antibodies can vary from 0.01% to 1%.

Our recommendation for usage:

- Use 150 µl **BSA Depletion Magnetic Agarose** (slurry) for 100 µl of antibody supplemented with $\leq 0.1\%$ BSA.
- Use 300 µl beads (slurry) for 100 µl antibody solutions with 0.5% BSA.
- Use 500 µl beads (slurry) for 100 µl antibody solutions with 1% BSA.

Bead volumes can be scaled up proportionally to the antibody volume.

- Please note: Commercial BSA preparations often contain bovine serum proteins, including bovine immunoglobulins (IgG). As a result, the recovered IgG may appear lower in quantity, when in fact the purification has removed contaminating bovine IgG and yielded a cleaner preparation of the target IgG.

Procedure

Equilibration of beads (for 100 µl antibody volume)

1. Resuspend beads by pipetting up and down or by inverting. **Do not vortex.**
2. Transfer 200 µl (for BSA concentrations < 0.5%) or 500 µl (for BSA 1%) beads (slurry) to a 1.5 mL tube. Separate beads with a magnet (or centrifuge at 2000x g for 1 min). Discard supernatant.
3. Add 500 µl PBS to the beads and separate beads with a magnet (or centrifuge at 2000x g for 1 min). Discard supernatant.

Depletion of BSA

4. Add your antibody (100 µl) to the equilibrated beads.
5. Incubate for 15 min at room temperature at 1000 rpm in a thermal shaker.
6. Separate beads with a magnet (or centrifuge at 2000x g for 1 min). The supernatant contains your BSA depleted antibody.
7. Optional (may increase antibody yield): Add 100 µl PBS to the beads and Separate beads with a magnet (or centrifuge at 2000x g for 2 min). Combine with supernatant from step 6.

Note: This will decrease the concentration of the antibody.