# **RFP-Trap for Immunoprecipitation**

From cell to gel in less than one hour

- ► No heavy & light chains in your downstream application
- Robust performance with high yield
- Short incubation time of about 30 minutes



# RFP-Trap®

ChromoTek RFP-Trap is a fully validated, well-established tool for fast and efficient one-step immunoprecipitation of RFP-fusion protein variants and their interacting factors from cell extracts or organelles. For research applications only.

# **Applications**


- Immunoprecipitation (IP) / Co-IP
- Mass spectrometry
- Enzyme activity measurements
- ChIP / RIP analysis

## Recognizes common RFP variants

- Including mCherry, mKate2, mPlum, mRFP, mRFPruby, tagRFP etc.
- Complete list at www.chromotek.com

# **Technology**

Camelidae, e.g. alpacas, possess a type of antibodies called heavy chain antibodies. These are devoid of light chains and bind their antigen via a single variable domain ( $V_HH$ ), also known as a nanobody. RFP-Trap consists of an alpaca single domain antibody that is coupled to agarose, magnetic agarose, and magnetic beads.



RFP-Trap: anti-RFP  $V_HH$ coupled to beads


# **Properties**

- Highly affine with a low dissociation constant of 4nM
- Very high binding capacity of more than 8µg RFP per 10µl slurry; magnetic beads less
- Stable in harsh conditions up to 2M NaCl or 2M urea
- Versatile: coupled to agarose, magnetic agarose, and magnetic beads

# Validated affinity reagent

- Structure and function are fully characterized
- Produced at constant high quality virtually without batch-to-batch variation
- Proven: used for experiments published in more than
  50 scientific articles

## Single band purification



Effective immunoprecipitations from protein extracts of human cells. Input (I), non-bound (FT) and bound fractions of mRFP, mCherry & mOrange were separated by SDS-PAGE followed by Coomassie staining.

# anti-RFP IgG Antibody 5F8 for IF

The most frequently cited anti-RFP antibody

- First in class: more than 100 publications
- Validated rat monoclonal antibody subtype IgG2a

# anti-mRFP IgG antibody 5F8

ChromoTek anti-RFP antibody 5F8 is a reliable tool for sensitive detection of RFP-fusion protein variants in immunofluorescence (IF). It is ranked #1 referenced anti-RFP antibody at www.citeAB.com. For research applications only.

IF of HeLa cells stably expressing mCherrys-PCNA

# Specificity

Tested on dsRed, mCherry, mPlum, mRFP, mRFPruby

## **Applications**

Immunofluorescence, ELISA, Immunoprecipitation


# anti-RFP IgG Antibody 6G6 for WB

**Exceptional detection of red fluorescent protein** variants

- Very sensitive: down to 4ng RFP / band
- Validated mouse monoclonal antibody subtype IgG2c

# anti-RFP IgG antibody 6G6

ChromoTek anti-RFP antibody 6G6 is an exceptional tool for sensitive detection of RFP-fusion protein derivatives in Western blots (WB). For research applications only.



WB c trans red fl

# **Specificity**

Tested on dsRed, mCherry, mOrange, mPlum, mRFP, mRFPruby

### **Applications**

Western blot

### ChromoTek Inc., USA

400 Oser Ave Ste 1650

Hauppauge, NY 11788

Phone: 631. 501. 1058 Fax: 631. 501. 1060

All other orders:

**Order information** 

coupled to agarose beads

Size

10 rxns

20 rxns

100 rxns

200 rxns

400 rxns

20 rxns

10 rxns

20 rxns

100 rxns

200 rxns

400 rxns

20 rxns

10 rxns

20 rxns

100 rxns

200 rxns

400 rxns

20 rxns

10 units

20 units

50 units

20 rxns

20 rxns

20 rxns

250µg

20µl

100µl

20µl

100µl

Order online - Save shipping costs!

USAorders@chromotek.com ➤ fax 631.501.1060

order@chromotek.com ▶ fax +49. 89.124.148.811

Code

rta-10

rta-20

rta-100

rta-200

rta-400

rtak-20

rtma-10

rtma-20

rtma-100

rtma-200

rtma-400

rtmak-20

rtm-10

rtm 20

rtm-100

rtm-200

rtm-400

rtmk-20

sct-10

sct-20

sct-50

bab-20

bmab-20

bmp-20

rt-250

5F8-20

5F8-100

6G6-20

6G6-100

Product name &

description

RFP-Trap A

▶ anti-RFP V.H.

RFP-Trap A Kit

▶ incl. lysis, wash and elution buffers

coupled to magnetic

► RFP-Trap A

RFP-Trap MA

agarose beads

RFP-Trap MA Kit

incl. lysis, wash and elution buffers

coupled to magnetic

► RFP-Trap MA

RFP-Trap M

► anti-RFP V<sub>H</sub>H

agarose beads

RFP-Trap M Kit

incl. lysis, wash and elution buffers

► RFP-Trap M

Spin columns

Binding control

Binding control

Binding control

▶ magnetic agarose beads

► magnetic agarose beads anti-RFP antibodies:

► Camelid; purified V\_H protein for IP, ELISA, and IF.

anti-RFP antibody [5F8]

anti-RFP antibody [6G6]

online www.chromotek.com/webshop

online www.chromotek.com/webshop

Mouse monoclonal IgG2c

► Rat monoclonal IgG2a for

Immunofluorescence

for Western Blotting

Used in REP-Trap

anti-RFP V<sub>H</sub>H/RFP binding protein

▶ agarose beads

► anti-RFP V<sub>H</sub>H

ell extract of HEK cells	
ently expressing different	
uorescent proteins	

mPlum- mRFP-PCNA PCNA morange mCherry dsRed TagRFP

# ChromoTek GmbH, Germany

Am Klopferspitz 19

D-82152 Planegg-Martinsried

Phone: +49. 89. 124. 148. 80 Fax: +49. 89. 124. 148. 811