

# Spin column protocol for Immunoprecipitation and Glycine elution of proteins from Nano-Trap<sup>®</sup> agarose beads

Only for research applications, not for diagnostic or therapeutic use.

<b>Harvest cells</b>	<p>For one immunoprecipitation reaction the use of <math>\sim 10^6</math> - <math>10^7</math> mammalian cells (approx. one 10-cm dish) expressing a protein of interest is recommended. To harvest adherent cells, aspirate growth medium, add 1 ml ice-cold PBS to cells and scrape cells from dish. Transfer cells to a pre-cooled tube, spin at 500 g for 3 min at +4°C and discard supernatant. Wash cell pellet twice with ice-cold PBS, gently resuspending the cells. After washing:</p>
<b>Lyse cells</b>	<ol style="list-style-type: none"><li>1. Resuspend cell pellet in 200 <math>\mu</math>l ice-cold lysis buffer by pipetting or using a syringe. <i>note: Supplement lysis buffer with protease inhibitors and 1 mM PMSF (not included).</i> <i>optional for nuclear/chromatin proteins: Use RIPA buffer supplemented with 1 mg/ml DNase, 2.5 mM MgCl<sub>2</sub>, protease inhibitors and 1 mM PMSF (not included).</i></li><li>2. Place the tube on ice for 30 min with extensively pipetting every 10 min.</li><li>3. Centrifuge cell lysate at 20.000x g for 10 min at +4°C. Transfer supernatant to a pre-cooled tube. Add 300 <math>\mu</math>l dilution buffer to lysate. Discard pellet. <i>note: At this point cell lysate may be put at -80°C for long-term storage.</i> <i>optional: Add 1 mM PMSF and protease inhibitors (not included) to dilution buffer.</i></li></ol>
<b>Equilibrate beads</b>	<ol style="list-style-type: none"><li>4. Remove the upper cap of a new spin column and snap of the plug from the bottom of the spin column (Keep cap and plug!). Place the spin column in a 2 ml tube. Vortex Nano-Trap<sup>®</sup>_A beads and pipette 25 <math>\mu</math>l bead slurry into the spin column. Immediately add 500 <math>\mu</math>l ice-cold dilution buffer. Centrifuge at 100x g for 5-10 sec. Discard flow-through and repeat wash twice. Close column with the bottom plug.</li></ol>
<b>Bind proteins</b>	<ol style="list-style-type: none"><li>5. Add diluted lysate (step 3) to equilibrated Nano-Trap<sup>®</sup>_A beads (step 4). If required, save 50 <math>\mu</math>l of diluted lysate for immunoblot analysis. Secure the top of the spin column. Tumble end-over-end for 1 hour at 4°C.</li><li>6. Remove the bottom cap from the spin column and place it in a new 2 ml tube. Centrifuge at 100x g for 5-10 sec. If required, save 50 <math>\mu</math>l flow-through for immunoblot analysis. Discard remaining flow-through.</li></ol>
<b>Wash beads</b>	<ol style="list-style-type: none"><li>7. Resuspend Nano-Trap<sup>®</sup>_A beads in 500 <math>\mu</math>l ice-cold dilution buffer. Place spin column in a 2 ml tube and centrifuge at 100x g for 5-10 sec. Discard flow-through and repeat wash twice. Close column with the bottom plug. <i>optional: Increase salt concentration in the second washing step up to 500 mM.</i></li></ol>
<b>Elute proteins</b>	<ol style="list-style-type: none"><li>8. Add 50 <math>\mu</math>l elution buffer (0.2 M glycine, pH 2.5) to Nano-Trap<sup>®</sup>_A beads. Pipette Nano-Trap<sup>®</sup>_A beads up and down for 30 sec. Remove bottom plug of the spin column and place it in a new 2 ml tube containing 5 <math>\mu</math>l 1M Tris base pH 10.4. Centrifuge at 1000x g for 30-60 sec. To increase elution efficiency this step can be repeated.</li></ol>

**Support** Please refer to our FAQ section at [www.chromotek.com](http://www.chromotek.com) or contact [support@chromotek.com](mailto:support@chromotek.com)