

NEURONAL MARKER ANTIBODIES

www.ptglab.com

Focus On WFS1

Catalog Number
11558-1-AP

Applications
CoIP, ELISA, IF, IHC, IP, WB

18 Publications

si Tested with siRNA

Wolfram syndrome protein (WFS1), also known as Wolframin, is a ubiquitously expressed protein with the highest levels of expression found in the brain, pancreas, heart, and insulinoma beta-cell lines. WFS1 is involved in the regulation of calcium homeostasis, which aids in the maintenance of the endoplasmic reticulum. Mutations in the WFS1 gene are associated with Wolfram syndrome, an autosomal recessive neurodegenerative disorder characterized by diabetes mellitus, optic atrophy, and various additional problems involving the urinary tract, brain, and hearing.

Island and Ocean Cells

Neuroscientists based at MIT (Massachusetts Institute of Technology), led by Nobel Prize-winning scientist Susumu Tonegawa, have identified and characterized two types of


cell in the entorhinal cortex that process the 'when' and 'where' signals in the brain - 'island' and 'ocean' cells.

Proteintech's WFS1 Antibody in Action

Island and ocean cells are distinguishable by their expression of certain neuronal markers. Ocean cells are identifiable as reelin-positive stellate cells, (Fig.B), whereas island cells are pyramidal neurons that stain positive for WFS1, (Fig.A). Though wolframin expression is

ubiquitous throughout the body, Tonegawa's group found that it is only present in island cells in the entorhinal cortex. The group used Proteintech's anti-WFS1 antibody (11558-1-AP) to selectively stain island cells within their research, as shown in several papers.

Download this catalog and lots like it at www.ptglab.com.



Golgi stained examples of A) a pyramidal cell and B) a stellate cell. (Churchill et al. *BMC Neuroscience* 2004 5:43 doi:10.1186/1471-2202-5-43 Figure 1).

Staining the two populations of cells with their specific markers allowed the Tonegawa group to map the trajectory of their connections to other parts of the brain. It became clear that island cells and ocean cells, though close in

proximity, had different targets: island cells plug into the CA1 region of the hippocampus, whereas ocean cells connect to the dentate gyrus and CA3 region of the hippocampus.

Focus On IBA1

Catalog Number
10904-1-AP

Applications
ELISA, IF, IHC, IP, WB

38 Publications

si Tested with siRNA

Ionized calcium binding adapter molecule 1 (IBA1), is an inflammatory response protein expressed by monocytes and macrophages and is primarily considered a marker of macrophage activation.

IBA1 is commonly used as a marker for microglial cells. Microglial cells are resident macrophages found in the CNS and play an important role in modulating inflammatory responses within the brain. Microglia are known to scavenge both the brain and spinal cord for pathogens and cells undergoing apoptosis and aid in their removal through phagocytosis.

Activation of microglia are one of the characteristic cellular changes that occurs during neurodegenerative processes such as Parkinson's disease and Alzheimer's disease. A recent study by Cui et al. has examined the etiological association between noise exposure and Alzheimer's disease. It has long been known that Alzheimer's disease results in neuroinflammation, therefore, Cui et al. examined the level of neuroinflammation present in the brains of mice exposed to chronic noise. Using Proteintech's antibody IBA1 as a marker for microglia, the authors were able to detect neuroglial activation during neuroinflammation through Western blot analysis.

Furthermore, Cui et al. were able to demonstrate that when the mice were


exposed to noise, the expression of IBA1 significantly increased for up to 14 days suggesting that chronic noise exposure may induce profound microglia activation, which may play a role in neuroinflammation and neurodegeneration.

As well as playing an important role in neurodegenerative disorders, IBA1 is also thought to be involved in neuropathic pain, a condition characterised by damage to the somatosensory nervous system.


A recent study by Popiolek-Barczyk et al. examined the use of the compound parthenolide on the ability to activate M2 microglia/macrophages following sciatic nerve injury. The authors demonstrated that parthenolide increased the levels of IBA1 using our antibody alongside a reduction in pain.

IBA1 clearly plays an important role in a variety of neurological conditions including neurodegenerative processes and pain. The Proteintech IBA1 antibody has been extensively validated in-house, including specificity testing using siRNA. Furthermore, it has been used in 25 peer-reviewed publications.

You can place your order online, by email or fax.



Immunohistochemical of paraffin-embedded human lymphoma using IBA1 antibody (10904-1-AP) at dilution of 1:100 (under 40x lens).



Immunofluorescent analysis of fixed C6 cells (-20°C Ethanol) using 10904-1-AP (IBA1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG (H+L).

References

- Cui B, Li K, Gai Z, et al. Chronic Noise Exposure Acts Cumulatively to Exacerbate Alzheimer's Disease-Like Amyloid- β Pathology and Neuroinflammation in the Rat Hippocampus. *Scientific Reports*. 2015;5:12943. doi:10.1038/srep12943.
- Popiolek-Barczyk K, Kolosowska N, Piotrowska A, et al. Parthenolide Relieves Pain and Promotes M2 Microglia/Macrophage Polarization in Rat Model of Neuropathy. *Neural Plasticity*. 2015;2015:676473. doi:10.1155/2015/676473.

THE BENCHMARK IN VALIDATION

Proteintech take complete responsibility for every product they manufacture and ensure that it will work in your hands with an extensive validation process.

Here's what you can expect:

- **Every antibody has been extensively validated by the in-house team in a wide range of applications.**
- **Antibodies now tested in siRNA knockdown to demonstrate specificity.**
- **Original validation data available to researcher on request.**
- **Multi-lingual technical support available online and by email.**

www.ptglab.com



ABAT
→ Galectin-3

ANTIBODY PRODUCT LIST

Antibody Name	Cat. No.	Type	Applications
ABAT	1 11349-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
ABCA2	20681-1-AP	Rabbit Poly	ELISA, IF, WB
ADH1B	1 17165-1-AP	Rabbit Poly	ELISA, WB
ALDH1L1	17390-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
Aquaporin 4	3 16473-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
BCL11B	55414-1-AP	Rabbit Poly	ELISA, IF, IP, WB
BCRP/ABCG2	2 10051-1-AP	Rabbit Poly	ELISA, IHC, WB
beta Tubulin	24 10094-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
beta Tubulin	17 10068-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
beta Tubulin	2 66240-1-Ig	Mouse Mono	ELISA, IF, IHC, IP, WB
BRN2	1 14596-1-AP	Rabbit Poly	ELISA, IHC, WB
BRN2-Specific	18998-1-AP	Rabbit Poly	ELISA, IHC, WB
Calretinin	12278-1-AP	Rabbit Poly	ELISA, IHC, WB
CAMSAP1L1	si 9 17880-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
CD11B/ Integrin alpha M	2 20991-1-AP	Rabbit Poly	ELISA, FC, IHC, WB
CD11B/ Integrin alpha M	21851-1-AP	Rabbit Poly	ELISA, IHC
CD133	18495-1-AP	Rabbit Poly	ELISA, IHC, WB
CD133	si 17 18470-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, WB
CD133-1, 2	19946-1-AP	Rabbit Poly	ELISA, IP, WB
CD133-1, 2, 4, 7	19945-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
CD24	2 18330-1-AP	Rabbit Poly	ELISA, WB
CD24	10600-1-AP	Rabbit Poly	ELISA, IHC, WB
CHAT	3 20747-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
CHAT	24418-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
CHT1	21848-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
Claudin 11	1 12152-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, WB


00 This number shows the amount of times our antibody has been cited in a publication.

Antibody Name	Cat. No.	Type	Applications
CNPase	1 13427-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
CUX1	si 2 11733-1-AP	Rabbit Poly	ELISA, IF, IP, WB
CUX2	24902-1-AP	Rabbit Poly	ELISA, WB
DACH1	si 26 10914-1-AP	Rabbit Poly	ChIP, ELISA, IF, IHC, IP, WB
DAT	3 22524-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
DCX	13925-1-AP	Rabbit Poly	ELISA, WB
DOPA decarboxylase	2 10166-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
DRD1	si 2 17934-1-AP	Rabbit Poly	ELISA, WB
DRD2	si 1 55084-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
DRD5	si 3 20310-1-AP	Rabbit Poly	ELISA, FC, WB
DYNC1H1	si 9 12345-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
EAAT2	22515-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
EAAT3	12686-1-AP	Rabbit Poly	ELISA, IHC, WB
EAAT4	12876-1-AP	Rabbit Poly	ELISA, WB
FABP7	17456-1-AP	Rabbit Poly	ELISA, WB
FABP7	1 51010-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
FABP7-Specific	14836-1-AP	Rabbit Poly	ELISA, IHC, WB
Ferritin light chain	si 7 10727-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
FEV	25058-1-AP	Rabbit Poly	ELISA, WB
FOXA1	20411-1-AP	Rabbit Poly	ELISA, IHC, WB
FOXA2	22474-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
FOXP2	20529-1-AP	Rabbit Poly	ELISA, IP, WB
GAC-specific	6 19958-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
GAD2	21760-1-AP	Rabbit Poly	ELISA, IP, WB
GAD2	1 20746-1-AP	Rabbit Poly	ELISA, IF, IP, WB
Galectin-3	si 2 14979-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB

si This icon shows the antibody has been tested in siRNA treated samples.

Galectin-3
PAX3 ←

More validation images available on our website.



Antibody Name	Cat. No.	Type	Applications
Galectin-3	1 60207-1-Ig	Mouse Mono	ELISA, IF, IHC, WB
GATA2	2 11103-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
GATA3	8 10417-1-AP	Rabbit Poly	ChIP, ELISA, IHC, WB
GATA3	22343-1-AP	Rabbit Poly	ELISA, WB
GFAP	23935-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
GFAP	22 16825-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
GFAP	6 60190-1-Ig	Mouse Mono	ELISA, IF, IHC, IP, WB
GIRK2	21647-1-AP	Rabbit Poly	ELISA, WB
GLAST	20785-1-AP	Rabbit Poly	ELISA, IHC, WB
Glutamine synthetase	4 11037-2-AP	Rabbit Poly	ELISA, FC, IHC, WB
GOT1	3 14886-1-AP	Rabbit Poly	ELISA, IF, IP, WB
GOT1	60317-1-Ig	Mouse Mono	ELISA, IHC, WB
HES5	22666-1-AP	Rabbit Poly	ELISA, IHC, WB
HMGB1	24 10829-1-AP	Rabbit Poly	ChIP, ELISA, FC, IF, IHC, IP, WB
HTR3A	1 10443-1-AP	Rabbit Poly	ELISA, FC, IHC, WB
HuC/D	13032-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
HuD	24992-1-AP	Rabbit Poly	ELISA, IF, IP, WB
IBA1	38 10904-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
Integrin beta-4	21738-1-AP	Rabbit Poly	ELISA, FC, IHC, IP, WB
Islet 1	15661-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
KGA/GAC	23 12855-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB Wes System
KGA/GAM/GAC	1 23549-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
KGA-Specific	3 20170-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
KIF5B	1 21632-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
LHX5	21567-1-AP	Rabbit Poly	ELISA, WB
LMX1B	18278-1-AP	Rabbit Poly	ELISA, IF, WB
MAG	14386-1-AP	Rabbit Poly	ELISA, IHC, WB
MAP1B	3 21633-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, WB
MAP2	17 17490-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
MOG	3 12690-1-AP	Rabbit Poly	ELISA, Inhibition Assay, WB


Antibody Name	Cat. No.	Type	Applications
Myelin basic protein	2 10458-1-AP	Rabbit Poly	ELISA, IHC, WB
NANOG	5 14295-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
Nestin	7 19483-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
NeuN	1 23060-1-AP	Rabbit Poly	ELISA, IHC, WB
NF200	60331-1-Ig	Mouse Mono	ELISA, IF, IHC, WB
NF-H	6 18934-1-AP	Rabbit Poly	ELISA, IHC, WB
NF-H	21471-1-AP	Rabbit Poly	ELISA, IHC, WB
NFIA	11750-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
NF-L	60189-1-Ig	Mouse Mono	ELISA, IHC, WB
NF-L	1 12998-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
NF-M-Specific	20664-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
NKX2-2	3 13013-1-AP	Rabbit Poly	ELISA, IF, WB
NOTCH1	3 20687-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
Notch1	5 10062-2-AP	Rabbit Poly	ELISA, IHC, IP, WB
NR2B	2 19954-1-AP	Rabbit Poly	ELISA, FC, IF, WB
NR2B	4 21920-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
NR2F1	24573-1-AP	Rabbit Poly	ELISA, IHC, WB
NSE	55235-1-AP	Rabbit Poly	ELISA, IHC, WB
NSE	10149-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
NSE	66150-1-Ig	Mouse Mono	ELISA, IF, IHC, WB
Nurr1/NR4A2	2 10975-2-AP	Rabbit Poly	ELISA, IHC, WB
OCT4	60242-1-Ig	Mouse Mono	ELISA, IHC, WB
OCT4	20 11263-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
Oct6-Specific	18997-1-AP	Rabbit Poly	ELISA, WB
OLIG2	4 13999-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
OLIG2	25754-1-AP	Rabbit Poly	ELISA, WB
p75NTR	2 55014-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
PAX2	2 21385-1-AP	Rabbit Poly	ELISA, IF, WB
PAX3	21386-1-AP	Rabbit Poly	ELISA, IP, WB
PAX3	51036-2-AP	Rabbit Poly	ELISA, WB









00 This number shows the amount of times our antibody has been cited in a publication.

si This icon shows the antibody has been tested in siRNA treated samples.











PAX7
→ WFS1


More validation images available on our website. 



Antibody Name	Cat. No.	Type	Applications
PAX7	20570-1-AP	Rabbit poly	ELISA, WB
Peripherin	17399-1-AP	Rabbit poly	ELISA, IHC, WB
PHOX2B  	25276-1-AP	Rabbit poly	ELISA, IF, IHC, WB
PHOX2B	66254-1-Ig	Mouse Mono	ELISA, IHC, WB
PLAGL1	25864-1-AP	Rabbit poly	ELISA, WB
RELN	20689-1-AP	Rabbit poly	ELISA, IHC
RORB	17635-1-AP	Rabbit poly	ELISA, IHC, IP, WB
S100B  	15146-1-AP	Rabbit poly	ELISA, IF, IHC, WB
SCN10A 	55334-1-AP	Rabbit poly	ELISA, IF
SCN9A/ Nav1.7-Specific 	20257-1-AP	Rabbit poly	ELISA, IF, IHC, WB
Somatostatin -Specific	17512-1-AP	Rabbit poly	ELISA, IHC, WB
SOX10	10422-1-AP	Rabbit poly	ELISA, WB
SOX2  	20118-1-AP	Rabbit poly	ELISA, IHC, WB

 This number shows the amount of times our antibody has been cited in a publication.

Antibody Name	Cat. No.	Type	Applications
SOX2 	11064-1-AP	Rabbit poly	ELISA, IF, IHC, WB
TAC1	13839-1-AP	Rabbit poly	ELISA, IHC, WB
TBR1 	20932-1-AP	Rabbit poly	ELISA, IF, IHC, IP, WB
TH 	25859-1-AP	Rabbit poly	ELISA, IF, IHC, IP, WB
TPH1	20879-1-AP	Rabbit poly	ELISA, WB
TPH2	22590-1-AP	Rabbit poly	ELISA, WB
TrkB	13129-1-AP	Rabbit poly	ELISA, IHC, IP, WB
VGLUT1	55491-1-AP	Rabbit poly	ELISA, FC, IHC
VGLUT2 	25261-1-AP	Rabbit poly	ELISA, IF, WB
Vimentin	22031-1-AP	Rabbit poly	ELISA, IHC, WB
Vimentin  	60330-1-Ig	Mouse Mono	ELISA, IF, IHC, WB
Vimentin  	10366-1-AP	Rabbit poly	ELISA, FC, IF, IHC, WB
WFS1  	11558-1-AP	Rabbit poly	ColP, ELISA, IF, IHC, IP, WB

 This icon shows the antibody has been tested in siRNA treated samples.



THE BENCHMARK IN ANTIBODIES

Proteintech products now feature in over
20,000 publications worldwide.

www.ptglab.com



CONTACT US

Proteintech Group

US Head Office

PHONE 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

FAX 1 (312) 455-8408

ADDRESS Proteintech Group, Inc.
5400 Pearl Street, Suite 300,
Rosemont, IL 60018, USA

EMAIL proteintech@ptglab.com

Proteintech Europe

United Kingdom

PHONE +44 (161) 8393007

FAX +44 (161) 2413103

ADDRESS Proteintech Europe, Ltd.
4th Floor,
196 Deansgate,
Manchester, M3 3WF

EMAIL europe@ptglab.com

Proteintech Europe

Germany

EMAIL germany@ptglab.com
Sales and technical support only.

Proteintech

China Office

PHONE 027-87531629
or 4006-900-926

FAX 027-87531627

ADDRESS Wuhan Sanying Biotechnologies
D3-3, No.666 Gaoxin Avenue,
Wuhan East Lake Hi-tech Development Zone
Wuhan, Hubei, P.R.C

EMAIL service@ptglab.com

Support

Available 24 hours via Live Chat and 9-5 (CDT) via phone.

Email support also available.

LIVE CHAT
www.ptglab.com

TWITTER
[@proteintech](https://twitter.com/proteintech)

BLOG
www.ptglab.com/news/blog

YOUTUBE
www.youtube.com/Proteintech



We are ISO 9001 and
ISO 13485 accredited.

