For Research Use Only

## CoraLite® Plus 647-conjugated BCL2 Recombinant antibody

proteintech®
Antibodies | ELISA kits | Proteins
www.ptglab.com

Catalog Number: CL647-80313

**Basic Information** 

Catalog Number:

CL647-80313

Size: 1000 µg/ml Source: Rabbit

Isotype: IgG

Immunogen Catalog Number:

AG3508

GenBank Accession Number: BC027258

GeneID (NCBI):

UNIPROT ID: P10415 Full Name:

B-cell CLL/lymphoma 2

Calculated MW: 26 kDa

Observed MW: 26 kDa

Purification Method: Protein A purification

CloneNo.:

305

Recommended Dilutions:

IF 1:50-1:500

Excitation/Emission maxima

wavelengths: 654 nm / 674 nm

**Applications** 

Tested Applications: FC (Intra), IF-P Species Specificity:

Human

Positive Controls:

IF: human appendicitis tissue,

## **Background Information**

BCL2 belongs to the Bcl-2 family. It suppresses apoptosis in a variety of cell systems including factor-dependent lymphohematopoietic and neural cells. BCL2 regulates cell death by controlling the mitochondrial membrane permeability (PMID: 10365962). It appears to function in a feedback loop system with caspases. BCL2 inhibits caspase activity either by preventing the release of cytochrome c from the mitochondria and/or by binding to the apoptosis-activating factor (APAF-1) (PMID: 10444588).

Storage

Storage:


Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer


PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human appendicitis tissue using Coralite® Plus 647 BCL2 antibody (CL647-80313, Clone: 305) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human appendicitis tissue using Coralite® Plus 647 BCL2 antibody (CL647-80313, Clone: 305) at dilution of 1:200.



1X10^6 Jurkat cells were intracellularly stained with 0.4 ug Coralite® Plus 647 Anti-Human BCL2 (CL647-80313, Clone:305) (red), or 0.4 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).