

For Research Use Only

HAS1 Polyclonal antibody

Catalog Number: 15229-1-AP **1 Publications**



Basic Information

Catalog Number: 15229-1-AP	GenBank Accession Number: BC035837	Purification Method: Antigen affinity purification
Size: 450 µg/ml	GeneID (NCBI): 3036	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q92839	
Isotype: IgG	Full Name: hyaluronan synthase 1	
Immunogen Catalog Number: AG4901	Calculated MW: 65 kDa	
	Observed MW: 65 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse testis tissue,
Cited Applications: WB	
Species Specificity: human, mouse, rat	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xueqing Shen	30452312	Biosci Biotechnol Biochem	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 15229-1-AP (HAS1 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.