

For Research Use Only

NFKBIL1 Polyclonal antibody

Catalog Number: 26273-1-AP



Basic Information

Catalog Number:

26273-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21978

GenBank Accession Number:

BC105064

GeneID (NCBI):

4795

UNIPROT ID:

Q9UBC1

Full Name:

nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor-like 1

Calculated MW:

381 aa, 43 kDa

Observed MW:

40-43 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB: HepG2 cells, MDA-MB-453s cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 26273-1-AP (NFKB1L1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.