

For Research Use Only

# ZBTB24 Polyclonal ANTIBODY

Catalog Number: 27975-1-AP



## Basic Information

Catalog Number:  
27975-1-AP

Size:  
370 µg/ml

Source:  
Rabbit

Isotype:  
IgG

Purification Method:  
Antigen affinity purification

Immunogen Catalog Number:  
AC24458

GenBank Accession Number:  
BC036731

GeneID (NCBI):  
9841

Full Name:  
zinc finger and BTB domain containing 24

Calculated MW:  
697aa, 78 kDa

Observed MW:  
80-85 kDa

Recommended Dilutions:  
WB 1:500-1:3000

## Applications

Tested Applications:  
WB, ELISA

Species Specificity:  
Human

**Note: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Positive Controls:

WB : HeLa cells;

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
--------	-----------	---------	-------------

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27975-1-AP (ZBTB24 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours