

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant 419cn protein



---

## Basic Information

Catalog Number:

Ag24301

Size:

50 µg

Form:

Available lyophilized

Species:

Expression Source:

*e coli*-derived, PET30a, with N-terminal 6<sup>°</sup>His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

(419CN encoded by )

---

## Selected Validation Data

---

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.