


AKT SIGNALING PATHWAY

www.ptglab.com



Related Products

Legend	
—	Association
→	Directly Activates
----->	Indirectly Activates
—	Inhibits
<—	Binding
Ub	Ubiquitin
P	phosphorylates

抗体名称	货号	类型	应用	适用物种
AKT	60203-2-Ig	Mouse MAb	WB, IP, IHC, IF, ELISA	human, mouse, rat
AKT-phospho-S473	66444-1-Ig	Mouse MAb	WB, IHC, ELISA	human, mouse
AKT	10176-2-AP	Rabbit PAb	WB, IP, IHC, IF, FC, ELISA	human, mouse, rat, zebrafish
Caspase 9/p35/p10	66169-1-Ig	Mouse MAb	WB, IP, IHC, IF, ELISA	human, mouse
Caspase 9/p35/p10	10380-1-AP	Rabbit PAb	WB, IP, IHC, IF, FC, ELISA	human, mouse, rat
CREB1	12208-1-AP	Rabbit PAb	WB, IP, IHC, IF, FC, ELISA	human, mouse, rat, monkey
FKBP5	14155-1-AP	Rabbit PAb	WB, IP, IHC, ELISA	human, mouse, rat
FOXO1	66457-1-Ig	Mouse MAb	WB, ELISA	human
FOXO1	18592-1-AP	Rabbit PAb	WB, IHC, FC, ELISA	human, mouse, rat
FOXO3A	66428-1-Ig	Mouse MAb	WB, IP, IHC, IF, ELISA	human, mouse
FOXO3A	10849-1-AP	Rabbit PAb	WB, IHC, IF, ELISA	human, mouse, rat
FOXO4	19121-1-AP	Rabbit PAb	WB, IF, ELISA	human, mouse
FOXO4	21535-1-AP	Rabbit PAb	WB, IHC, ELISA	human, mouse, rat
GLUT4	66846-1-Ig	Mouse MAb	WB, ELISA	human, mouse, rat
GLUT4	21048-1-AP	Rabbit PAb	IHC, ELISA	human
GSK3A	13419-1-AP	Rabbit PAb	WB, IHC, IF, ELISA	human, mouse, rat, zebrafish
GSK3B	24198-1-AP	Rabbit PAb	WB, IP, IHC, IF, FC, ELISA	human, mouse
GSK3B	22104-1-AP	Rabbit PAb	WB, IP, IHC, IF, ELISA	human, mouse, rat

抗体名称	货号	类型	应用	适用物种
Hamartin	20988-1-AP	Rabbit PAb	WB, IP, IHC, IF, ELISA	human, mouse, rat
mTOR	66888-1-Ig	Mouse MAb	WB, IHC, ELISA	human, mouse, Rat
mTOR	20657-1-AP	Rabbit PAb	WB, IP, IHC, IF, ELISA	human, mouse, rat
p53	60283-2-Ig	Mouse MAb	WB, IHC, ELISA	human, mouse, rat
p53	10442-1-AP	Rabbit PAb	WB, IP, IHC, IF, CoIP, ChIP, ELISA	human, rat, mouse
p70(S6K)	66638-1-Ig	Mouse MAb	WB, IHC, ELISA	human, mouse, Rat
p70(S6K)	14485-1-AP	Rabbit PAb	WB, IP, IHC, IF, FC, ELISA	human, mouse, rat
Paxillin	10029-1-Ig	Rabbit PAb	WB, IP, IHC, IF, ELISA	human, mouse
PDK1	10026-1-AP	Rabbit PAb	WB, IP, IHC, ELISA	human, mouse
PI3K p110(beta)	20584-1-AP	Rabbit PAb	WB, IP, IHC, IF, co-IP, ELISA	human
PRAS40	21097-1-AP	Rabbit PAb	WB, IHC, IF, ELISA	human
PTEN	60300-1-Ig	Mouse MAb	WB, IHC, IF, ELISA	human, mouse
PTEN	22034-1-AP	Rabbit PAb	WB, IP, IHC, IF, ELISA	human, mouse
Raptor	20984-1-AP	Rabbit PAb	WB, IP, IHC, IF, ELISA	human, rat
RICTOR	27248-1-AP	Rabbit PAb	WB, IHC, IF, ELISA	human, mouse
RRAGC	26989-1-AP	Rabbit PAb	WB, ELISA	human, mouse
Tuberin	24601-1-AP	Rabbit PAb	WB, IHC, IF, ELISA	human
YBK1	20339-1-AP	Rabbit PAb	WB, IP, IHC, IF, ELISA	human, mouse

该数字表示抗体的文献引用次数