

THRA Polyclonal antibody

Catalog Number: 10139-1-AP

2 Publications

Basic Information

Catalog Number:

10139-1-AP

Size:

2500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0182

GenBank Accession Number:

BC002728

GeneID (NCBI):

7067

UNIPROT ID:

P10827

Full Name:

thyroid hormone receptor, alpha
(erythroblastic leukemia viral (v-erb-
a) oncogene homolog, avian)

Calculated MW:

54 aa, 7 kDa

Observed MW:

58 kDa, 50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

IHC, WB

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB : A431 cells, mouse brain tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Sijie Wang	36216077	Physiol Behav	WB
C B Steinhauser	34062290	Domest Anim Endocrinol	IHC

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 10139-1-AP (THRA antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.