

For Research Use Only

CARTPT Polyclonal antibody

Catalog Number: 13212-1-AP



Basic Information

Catalog Number:

13212-1-AP

Size:

187 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4004

GenBank Accession Number:

BC029882

GeneID (NCBI):

9607

UNIPROT ID:

Q16568

Full Name:

CART prepropeptide

Calculated MW:

116 aa, 13 kDa

Observed MW:

4-14 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : SH-SY5Y cells, human adrenal gland tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 13212-1-AP (CARTPT antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.