

For Research Use Only

# NOTUM Polyclonal antibody

Catalog Number: 14663-1-AP

2 Publications



## Basic Information

Catalog Number:

14663-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6338

GenBank Accession Number:

BC060882

GeneID (NCBI):

147111

UNIPROT ID:

Q6P988

Full Name:

notum pectinacetyltransferase homolog  
(Drosophila)

Calculated MW:

56 kDa

Observed MW:

56 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:10000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

rat

Positive Controls:

WB : mouse liver tissue,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yi-Feng Yuan	36304155	Front Pharmacol	WB
Panpan Yang	34352447	Arch Oral Biol	ELISA

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse liver lysates were subjected to SDS PAGE followed by western blot with 14663-1-AP (NOTUM antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.