

For Research Use Only

VPS37B Polyclonal antibody

Catalog Number: 15653-1-AP

Featured Product

5 Publications



Basic Information

Catalog Number:

15653-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8256

GenBank Accession Number:

BC005882

GeneID (NCBI):

79720

UNIPROT ID:

Q9H9H4

Full Name:

vacuolar protein sorting 37 homolog B
(S. cerevisiae)

Calculated MW:

31 kDa

Observed MW:

31 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human, mouse

Positive Controls:

WB : HeLa cells, COLO 320 cells

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|-------------------|-----------|------------|-------------|
| Simona M Migliano | 37840274 | Autophagy | WB |
| James Glover | 37439191 | J Cell Sci | WB |
| Ruheena Javed | 37272163 | EMBO J | WB |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 15653-1-AP (VPS37B Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.