

For Research Use Only

# ELF4 Polyclonal antibody

Catalog Number: 16020-1-AP **1 Publications**



## Basic Information

Catalog Number:

16020-1-AP

Size:

320 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8836

GenBank Accession Number:

BC017194

GeneID (NCBI):

2000

UNIPROT ID:

Q99607

Full Name:

E74-like factor 4 (ets domain transcription factor)

Calculated MW:

663 aa, 71 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Paul M Tyler	34326534	Nat Immunol	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data