

For Research Use Only

# Noxa Polyclonal antibody

Catalog Number: 17418-1-AP **2 Publications**



## Basic Information

Catalog Number:

17418-1-AP

Size:

270 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11387

GenBank Accession Number:

BC032663

GeneID (NCBI):

5366

UNIPROT ID:

Q13794

Full Name:

phorbol-12-myristate-13-acetate-induced protein 1

Calculated MW:

136 aa, 15 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Fangbing Liu	31115744	Target Oncol	WB
Huiying Zhao	28573788	Cancer Med	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data