

For Research Use Only

TSPAN32 Polyclonal antibody

Catalog Number: 20320-1-AP



Basic Information

Catalog Number: 20320-1-AP	GenBank Accession Number: BC016693	Purification Method: Antigen affinity purification
Size: 300 µg/ml	GeneID (NCBI): 10077	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	UNIPROT ID: Q96QS1	
Isotype: IgG	Full Name: tetraspanin 32	
Immunogen Catalog Number: AG14174	Calculated MW: 320 aa, 35 kDa	
	Observed MW: 33 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse spleen tissue,
Species Specificity: human, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse spleen tissue were subjected to SDS PAGE followed by western blot with 20320-1-AP (TSPAN32 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.