

For Research Use Only

TSPAN32 Polyclonal antibody

Catalog Number: 20320-1-AP



Basic Information

Catalog Number:

20320-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14174

GenBank Accession Number:

BC016693

GeneID (NCBI):

10077

UNIPROT ID:

Q96QS1

Full Name:

tetraspanin 32

Calculated MW:

320 aa, 35 kDa

Observed MW:

33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse spleen tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse spleen tissue were subjected to SDS PAGE followed by western blot with 20320-1-AP (TSPAN32 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.