

For Research Use Only

CRYGN Polyclonal antibody

Catalog Number: 21281-1-AP



Basic Information

Catalog Number: 21281-1-AP	GenBank Accession Number: BC100880	Purification Method: Antigen affinity purification
Size: 260 µg/ml	GeneID (NCBI): 155051	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	UNIPROT ID: Q8WXF5	
Isotype: IgG	Full Name: crystallin, gamma N	
Immunogen Catalog Number: AG15795	Calculated MW: 182 aa, 21 kDa	
	Observed MW: 17 kDa, 25 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse eye tissue,
Species Specificity: human, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse eye tissue were subjected to SDS PAGE followed by western blot with 21281-1-AP (CRYGN Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.