

For Research Use Only

BMP10 Polyclonal antibody

Catalog Number: 22858-1-AP



Basic Information

Catalog Number: 22858-1-AP	GenBank Accession Number: BC105063	Purification Method: Antigen Affinity purified
Size: 150 µg/ml	GeneID (NCBI): 27302	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: O95393	
Isotype: IgG	Full Name: bone morphogenetic protein 10	
Immunogen Catalog Number: AG18746	Calculated MW: 424 aa, 48 kDa	
	Observed MW: 48 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293 cells, MCF7 cells
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF7 cells were subjected to SDS PAGE followed by western blot with 22858-1-AP (BMP10 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



HEK-293 cells were subjected to SDS PAGE followed by western blot with 22858-1-AP (BMP10 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.