

For Research Use Only

BMP10 Polyclonal antibody

Catalog Number: 22858-1-AP



Basic Information

Catalog Number:

22858-1-AP

Size:

150 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18746

GenBank Accession Number:

BC105063

GeneID (NCBI):

27302

UNIPROT ID:

O95393

Full Name:

bone morphogenetic protein 10

Calculated MW:

424 aa, 48 kDa

Observed MW:

48 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, MCF7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF7 cells were subjected to SDS PAGE followed by western blot with 22858-1-AP (BMP10 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



HEK-293 cells were subjected to SDS PAGE followed by western blot with 22858-1-AP (BMP10 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.