

For Research Use Only

# ACD Polyclonal antibody

Catalog Number: 25849-1-AP **2 Publications**



## Basic Information

Catalog Number:

25849-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23021

GenBank Accession Number:

BC016904

GeneID (NCBI):

65057

UNIPROT ID:

Q96APO

Full Name:

adrenocortical dysplasia homolog  
(mouse)

Calculated MW:

58 kDa

Observed MW:

57-60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

IP, WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB: Jurkat cells,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Ji Hoon Lee	34321211	Sci Adv	IP
Hongmiao Hu	37624885	Sci Adv	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 25849-1-AP (ACD Antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.