

For Research Use Only

53BP2 Polyclonal antibody

Catalog Number: 26550-1-AP



Basic Information

Catalog Number:

26550-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24907

GenBank Accession Number:

BC040247

GeneID (NCBI):

7159

UNIPROT ID:

Q13625

Full Name:

tumor protein p53 binding protein, 2

Observed MW:

150 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:5000-1:20000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293T cells, MCF-7 cells, Jurkat cells, HeLa cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26550-1-AP (53BP2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.