

For Research Use Only

# EPB41L4A Polyclonal antibody

Catalog Number: 26596-1-AP



## Basic Information

Catalog Number:

26596-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23923

GenBank Accession Number:

BC031042

GeneID (NCBI):

64097

UNIPROT ID:

Q9HCS5

Full Name:

erythrocyte membrane protein band

4.1 like 4A

Observed MW:

36 kDa, 79 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:300-1:800

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : NIH/3T3 cells, Transfected HEK-293 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 26596-1-AP (EPB41L4A Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



WB result of EPB41L4A antibody (26596-1-AP; 1:300; room temperature for 1.5 hours) with negative control and EPB41L4A overexpressed Transfected HEK-293 cells.