

For Research Use Only

# MCT14 Polyclonal antibody

Catalog Number: 26953-1-AP



## Basic Information

Catalog Number:

26953-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25315

GenBank Accession Number:

BC065524

GeneID (NCBI):

151473

UNIPROT ID:

Q7RTX9

Full Name:

solute carrier family 16, member 14  
(monocarboxylic acid transporter 14)

Calculated MW:

510 aa, 56 kDa

Observed MW:

42-52 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, pig, mouse

Positive Controls:

WB : Caco-2 cells, pig kidney tissue, SH-SY5Y cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Caco-2 cells were subjected to SDS PAGE followed by western blot with 26953-1-AP (MCT14 antibody at dilution of 1:2000 incubated at room temperature for 1.5 hours.