

For Research Use Only

TANK Polyclonal antibody

Catalog Number: 27065-1-AP **3 Publications**



Basic Information

Catalog Number:

27065-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25740

GenBank Accession Number:

BC003388

GeneID (NCBI):

10010

UNIPROT ID:

Q92844

Full Name:

TRAF family member-associated

NFKB activator

Calculated MW:

14 kDa

Observed MW:

48 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human, rat

Positive Controls:

WB : THP-1 cells, HepG2 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Shumin Qin	28942222	Int Immunopharmacol	WB
Xuan Huang	31070748	J Mol Cell Biol	WB
Takako Okamoto	32580962	Cancer Res	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



THP-1 cells were subjected to SDS PAGE followed by western blot with 27065-1-AP (TANK Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.