

For Research Use Only

# PEX26 Polyclonal antibody

Catalog Number: 27472-1-AP

2 Publications



## Basic Information

Catalog Number:

27472-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26529

GenBank Accession Number:

BC016280

GeneID (NCBI):

55670

UNIPROT ID:

Q7Z412

Full Name:

peroxisomal biogenesis factor 26

Calculated MW:

34 kDa

Observed MW:

30-34 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : A431 cells, HEK-293 cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Richard G Lee	38129691	Nat Cell Biol	WB
Bianbian Yan	37957408	Dig Dis Sci	IHC, WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27472-1-AP (PEX26 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.