

For Research Use Only

FBXO11 Polyclonal antibody

Catalog Number: 27610-1-AP

Featured Product

2 Publications



Basic Information

Catalog Number:

27610-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24401

GenBank Accession Number:

BC012728

GeneID (NCBI):

80204

UNIPROT ID:

Q86XK2

Full Name:

F-box protein 11

Calculated MW:

927 aa, 104 kDa

Observed MW:

94-104 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, WB, ELISA

Cited Applications:

CoIP, WB

Species Specificity:

Human

Cited Species:

human

Positive Controls:

WB : LNCaP cells, U2OS cells

IP : LNCaP cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Lijiang Shao	32657545	FEBS Open Bio	WB
Hao Zhang	37837399	Clin Transl Med	WB, CoIP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27610-1-AP (FBXO11 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.



IP result of anti-FBXO11 (IP:27610-1-AP, 4ug; Detection:27610-1-AP 1:0) with LNCaP cells lysate 2400 ug.