

For Research Use Only

TGM7 Polyclonal antibody

Catalog Number: 27748-1-AP



Basic Information

Catalog Number:

27748-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26875

GenBank Accession Number:

BC153105

GeneID (NCBI):

116179

UNIPROT ID:

Q96PF1

Full Name:

transglutaminase 7

Observed MW:

70-75 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

Positive Controls:

WB : mouse liver tissue, NIH/3T3 cells, mouse colon tissue, rat liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27748-1-AP (TGM7 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.