

For Research Use Only

CDSN Polyclonal antibody

Catalog Number: 27953-1-AP **1 Publications**



Basic Information

Catalog Number: 27953-1-AP	GenBank Accession Number: BC031993	Purification Method: Antigen affinity purification
Size: 900 µg/ml	GeneID (NCBI): 1041	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	UNIPROT ID: Q15517	
Isotype: IgG	Full Name: corneodesmosin	
Immunogen Catalog Number: AG27688	Calculated MW: 528 aa, 52 kDa	
	Observed MW: 55-65 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse skin tissue,
Cited Applications: IF	
Species Specificity: Human, mouse	
Cited Species: mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yosuke Ishitsuka	35024686	JID Innov	IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 27953-1-AP (CDSN antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.