

For Research Use Only

FBXO30 Polyclonal antibody

Catalog Number: 28039-1-AP



Basic Information

Catalog Number:

28039-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27791

GenBank Accession Number:

BC024326

GeneID (NCBI):

84085

UNIPROT ID:

Q8TB52

Full Name:

F-box protein 30

Calculated MW:

745 aa, 82 kDa

Observed MW:

90 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HeLa cells, A549 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28039-1-AP (FBXO 30 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.