

For Research Use Only

CYP302a1 (mosquito) Polyclonal antibody



Catalog Number: 51097-1-AP

Basic Information

Catalog Number:

51097-1-AP

Size:

160 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0761

GenBank Accession Number:

AY947549

GeneID (NCBI):

Full Name:

Calculated MW:

59 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

mosquito

Background Information

CYP302a1, also named as disembodied (dib-Bm), belongs to the cytochrome P450 family. This antibody detects the CYP302a1 protein of mosquito.

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data