

For Research Use Only

# MUC4 Polyclonal antibody

Catalog Number: 55343-1-AP



## Basic Information

Catalog Number:

55343-1-AP

Size:

147 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_004532

GeneID (NCBI):

4585

UNIPROT ID:

Q99102

Full Name:

mucin 4, cell surface associated

Calculated MW:

232 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

MUC4, also named as HSA276359, ASGP, MUC-4, testis mucin and mucin-4 is a high molecular weight transmembrane mucin that is overexpressed in many carcinomas and is a risk factor associated with a poor prognosis. MUC4 have a key role in the progression of pancreatic cancer. MUC4 helps in the docking of tumor cells on the endothelial surface. This antibody is specific to MUC4.

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

