

For Research Use Only

# P53 Monoclonal antibody

Catalog Number: 60283-1-Ig

Featured Product

23 Publications



## Basic Information

Catalog Number:

60283-1-Ig

Size:

673  $\mu$ g/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0698

GenBank Accession Number:

BC003596

GeneID (NCBI):

7157

UNIPROT ID:

P04637

Full Name:

tumor protein p53

Calculated MW:

44 kDa

Purification Method:

Protein G purification

CloneNo.:

2D2F10

## Applications

Tested Applications:

ELISA

Cited Applications:

IF, IHC, IP, WB

Species Specificity:

human, mouse, rat

Cited Species:

human, rat, mouse

## Background Information

TP53, also named as P53 and NY-CO-13, belongs to the p53 family. It has 9 isoforms. In SDS-Page, the MW is about 53kd. TP53 acts as a tumor suppressor in many tumor types; induces growth arrest or apoptosis depending on the physiological circumstances and cell type. It is involved in cell cycle regulation as a trans-activator that acts to negatively regulate cell division by controlling a set of genes required for this process. The antibody reacts specifically with 53 kDa protein.

## Notable Publications

Author	Pubmed ID	Journal	Application
Fang Wu	34535139	Cancer Cell Int	WB
Kun Lu	29285195	Oncol Lett	IF
Xi Chen	33174023	Mol. Med Rep	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data