

For Research Use Only

# Anti-Mouse CD24 (M1/69)

Catalog Number: 65121-1-Ig



---

## Basic Information

Catalog Number:

65121-1-Ig

Size:

100ug, 500 µg/ml

Source:

Rat

Isotype:

IgG2b, kappa

GenBank Accession Number:

BC075622

GeneID (NCBI):

12484

Full Name:

CD24a antigen

Purification Method:

Affinity purification

CloneNo.:

M1/69

---

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

---

## Background Information

---

### Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide, pH 7.3.

---

For technical support and original validation data for this product please contact:

T: 4006900926


E: Proteintech-CN@ptglab.com

W: ptgcn.com

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> BALB/c mouse splenocytes were surface stained with 0.5 ug Anti-Mouse CD24 (65121-1-Ig, Clone: M1/69) and FITC anti-rat IgG2b Antibody at dilution 1:200 (green), or stained with 0.5 ug isotype control and FITC anti-rat IgG2b Antibody at dilution 1:200 (black). Cells were not fixed.