

For Research Use Only

MEA1 Monoclonal antibody

Catalog Number: 66472-1-Ig



Basic Information

Catalog Number: 66472-1-Ig	GenBank Accession Number: BC001754	Purification Method: Protein G purification
Size: 1500 µg/ml	GeneID (NCBI): 4201	CloneNo.: 1H8C2
Source: Mouse	Full Name: male-enhanced antigen 1	Recommended Dilutions: WB 1:500-1:2000
Isotype: IgG1	Calculated MW: 20 kDa	
Immunogen Catalog Number: AG0375	Observed MW: 28 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : rat testis tissue, human testis tissue
Species Specificity: human, rat	

Background Information

MEA1, also named as Male-enhanced antigen 1, is a 185 amino acid protein, which is highly expressed in testis. MEA1 may play an important role in spermatogenesis and/or testis development.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

