

For Research Use Only

# DTX2 Monoclonal antibody

Catalog Number: 67209-1-Ig



## Basic Information

Catalog Number:

67209-1-Ig

Size:

2400 µg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG12022

GenBank Accession Number:

BC064535

GeneID (NCBI):

113878

UNIPROT ID:

Q86UW9

Full Name:

deltex homolog 2 (Drosophila)

Calculated MW:

253aa, 28 kDa; 575aa, 64 kDa

Observed MW:

68 kDa

Purification Method:

Protein A purification

CloneNo.:

2D4F8

Recommended Dilutions:

WB 1:1000-1:6000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HeLa cells, HepG2 cells, MCF-7 cells, A549 cells

## Background Information

### Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67209-1-Ig (DTX2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.