

For Research Use Only

# ARL2BP Monoclonal antibody

Catalog Number: 67316-1-Ig



## Basic Information

Catalog Number:

67316-1-Ig

Size:

1600 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29261

GenBank Accession Number:

BC003087

GeneID (NCBI):

23568

UNIPROT ID:

Q9Y2Y0

Full Name:

ADP-ribosylation factor-like 2 binding protein

Calculated MW:

19 kDa

Observed MW:

21 kDa

Purification Method:

Protein G purification

CloneNo.:

1G5C2

Recommended Dilutions:

WB 1:1000-1:6000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB: A549 cells, HeLa cells, Jurkat cells, HEK-293 cells, HepG2 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67316-1-Ig (ARL2BP antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 67316-1-Ig (ARL2BP antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.