

For Research Use Only

AFP Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68686-1-PBS

Basic Information

Catalog Number:

68686-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC027881

GeneID (NCBI):

174

UNIPROT ID:

P02771

Full Name:

alpha-fetoprotein

Calculated MW:

69 kDa

Purification Method:

Protein G purification

CloneNo.:

1B1E9

Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP50072-1, AFP Monoclonal Matched Antibody Pair - BSA and Azide free measured by Sandwich ELISA. Capture antibody: 68686-1-PBS. Detection antibody: 68686-2-PBS (HRP conjugated). Standard: Eg0539. Range: 31.25-2000 pg/mL