

For Research Use Only

SALL4 Recombinant antibody, PBS Only (Detector)



Catalog Number:83039-3-PBS

Basic Information

Catalog Number: 83039-3-PBS	GenBank Accession Number: BC111714	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57167	CloneNo.: 230284D8
Source: Rabbit	UNIPROT ID: Q9UJQ4	
Isotype: IgG	Full Name: sal-like 4 (Drosophila)	
Immunogen Catalog Number: AG16076	Calculated MW: 1053 aa, 112 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00137-2, SALL4 Recombinant Matched Antibody Pair - BSA and Azide free. Capture antibody: 83039-2-PBS. Detection antibody: 83039-3-PBS. Standard: Ag16076. Range: 39-2500 pg/mL