

For Research Use Only

CA15-3,MUC1 Recombinant antibody, PBS Only (Detector)



Catalog Number:83311-1-PBS

Basic Information

Catalog Number: 83311-1-PBS	GenBank Accession Number: BC120975	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4582	CloneNo.: 240190A2
Source: Rabbit	UNIPROT ID: P15941	
Isotype: IgG	Full Name: mucin 1, cell surface associated	
Immunogen Catalog Number: AG20366	Calculated MW: 1264 aa, 123 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00335-1, CA15-3, MUC1 Recombinant Matched Antibody Pair - BSA and Azide Free. Capture antibody: 83311-3-PBS. Detection antibody: 83311-1-PBS. Standard: Ag20366. Range: 0.78-100 ng/mL.



Cytometric bead array standard curve of MP00335-2, CA15-3, MUC1 Recombinant Matched Antibody Pair - BSA and Azide Free. Capture antibody: 83311-2-PBS. Detection antibody: 83311-1-PBS. Standard: Ag20366. Range: 1.56-100 ng/mL.