

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

CDK5 Fusion Protein



Basic Information

Catalog Number:

Ag0686

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

DSCNGDLDP EIVKSF L FQ L L K G L G F C H S R N V L H R D L K P
Q N L L I N R N G E L K L A D F G L A R A F G I P V R C Y S A E V V T L W Y
R P P D V L F G A K L Y S T S I D M W S A G C I F A E L A N A G R P L F P
G N D V D D Q L K R I F R L L G T P T E E Q W P S M T K L P D Y K P Y P M
Y P A T T S L V N V V P K L N A T G R D L L Q N L L K C N P V Q
(92-273 aa encoded by BC005115)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.