

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

CYP302a1 Fusion Protein



Basic Information

Catalog Number:

Ag0761

Size:

50 µg

Form:

Available lyophilized

Species:

Expression Source:

e coli-derived, PET28a, with N-terminal 6^oHis.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MYLSPVKSIKISATLMTSRSCATVLENVKPNQIPGP
RGPFGLGNLYQYIPGIGKYSFDALHESGQDKYEKYGPI
VRETMVPGQDIVWLYDPNDIAAVLNDKTPGIYPSRRS
HTALAKYRRDRPNVYRTAGLLATNGIEWWKIRSELQK
GLSSPQSVRNFLPLTDKVTREFVASMNSTEHDCVPDF
MPAISRLNLELICVMAFDVRLDSFSDEQMKPNSLSSRL
MESAEVTNQSIPTDQGFQLWKYFETPAYRKLKRAQE
FMEKTAVELVSQKLYFDEDQKQLASGRHRSLSLEEY
LRNPNLELHDII GMAADLLLAGVHTSSYTTFALYHLCL
NPDAQDKLYQEA CRILPDPWECQIEAAALNSEASYC
RAVLKESLRLNPISIGVGRILNKDATLGGYHVPKGTVV
VTQNLVSCRQERYFKNPTKFI PERWMRETKEDEVNPLY
VLPFGHGMRSCIARRMAEQNMLVLLLRLIRSYEIDWKG
KVPMNIETKLINQPDQPIKIAFRSRKS
(1-515 aa encoded by AY947549)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.