

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

NCL Fusion Protein



Basic Information

Catalog Number:
Ag0859

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MVKLAKAGKNQGDPPKMAPPPKEVEEEDSEDEEMSED
EEDSSGEEVVIPIQKKGKAAATSAKKVVVSPTKKVA
VATPAKKA AVTPGKKAATPAKKT VTPAKAVTTPGK
KGATPGKALVATPGKGA AIPAGKAKNGKNAKKEDS
DEEEDDDSEDEEDEDDEDEDEIEPAAMKAAAAA
PASEDEDEDDEDDEDGKSKGIA YIEFKTEADAEKTF
EEKQGT EIDGRSISLYYTGEKGNQDYRGGKNSTWSG
ESKTLVLSNLSYSATEETLQEVFEKATFIKVPQNQNGK
SKGYAFIEFASFEDA KEALNSCNKREIEGRAIRLELQGP
RGSFNARSQPSKTLFVKGLSEDTTEETLKEFDG SVRA
RIVTDRETGSSKGFVDFNSEEDAKAAKEAMEDGEI
DGNKVTLDWAKPKGEGGFGRGGGRGGFGRGGGR
GGRGGFGGRGRGFGRGGFRGGGGGDHQPQ GK
KTKFE
(1-482 aa encoded by BC006494)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.