

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# JAK2 Fusion Protein



## Basic Information

Catalog Number:

Ag11264

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

EDRHQLPAPKWAELANLINNCMDYEPDFRPSFRAIIR  
DLNSLFTPDYELLTENDMLPNMRIGALGFSGAFEDRDP  
TQFEERHLKFLQQLGKGNFGSVEMCRYDPLQDNTGE  
VVAVKKLQHSTEEHLRDFEREIEILKSLQHDNIVKYKG  
VCYSAGRRNLKIMEYLPYGSLRDYQKHKERIDHIKLL  
QYTSQICKGMEYLGTKRYIHRDLATRNILVENENRVKI  
GDFGLTKVLPQDKEYYKVKEPGESPIFWYAPESLTESK  
FSVASDVWSFGVVLYELFTYIEKSKSPPAEFMRMIGN  
DKQGQMIVFHLLKNNGRRLPRPDGCPDEIYMIMTEC  
WNNNVNQRPFRDLALRVDQIRDNMAG  
(764-1129 aa encoded by BC039695)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.