

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

CLCN1 Fusion Protein



Basic Information

Catalog Number:

Ag13336

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

DSMILLGSVERSELQALLQRHLCPERRLRAAQEMARKL
SELPYDYGKARLAGEGLPGAPPGRPESEFAFVDEDEDED
LSGKSELPPSLALHPSTTAPLSPEEPNGPLPGHKQQPE
APEPAGQRPSIFQSLHCLLGRARPTKKITTDSTDLV
DNMSPEEIEAWEQEQLSQPVCFDSCCIDQSPFQLVE
QTTLHKHTLFLSLLGLHLAYVTSMGKLRGVLAEELQK
AIEGHTKSGVQLRPPLASFRNTTSTRKSTGAPSSAEN
WNLPEDRPGATGTGDVIAASPETVPSPSPPEPLSLAP
GKVEGELEELVESPGLEELADILQGPSLRSTDEEDE
DELIL

(644-988 aa encoded by BC112156)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.