

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

PRMT1 Fusion Protein



Basic Information

Catalog Number:

Ag1790

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MEVSCGQAESEKPAEDMTSKDYDFDSYAHFGIHE
EMLKDEVRTLTYRNSMFHNRHLFKDKVVLVDVGS GTI
LCMFAAKAGARKVIGIECSSISDYAVKIVKANKLDHV
VTIIKGGVEEVELPVEKV DIIISEWMGYCLFYESMLNTV
LYARDKWLPDGLIFPD RATLYVTAIEDRQYKDYKIHW
WENVYGFDMSCIKDVAIKEPLVDVVDPKQLVTNACLI
KEVDIYTVKVEDLFTSPFCLQV KRNDYVHALVA YFNI
EFTRCHKRTGFSTSPESPYTHWKQTVFYMEDYLVKT
GEEIFGTIGMRPNAKNNRDLDF TIDLDFKQQLCELSCS
TDYRMR

(1-343 aa encoded by BC019268)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.